



## Nr Kat. 4614

Zestaw z wiertłami powlekanymi prostymi



Kaseta jest lekka, kompaktowa z bezpiecznym zamknięciem. Wyposażona jest w silikonowe uchwyty zabezpieczające przed przemieszczaniem się elementów podczas transportu.



Możliwość sterylizacji w autoklawie (1000 procesów).

### Materiały:

Kaseta:  
żywica polyfenylosulfonowa Radel R – PPSU  
Rynienka - stal nierdzewna

### Wymiary kasety:


Kaseta duża: 19 cm x 14cm x 6cm  
Nr kat. 4613



4 wygodne pozycje ustawienia. Laserowo wykonane oznaczenia są trwałe i wygodne do czyszczenia.



Kaseta chirurgiczna zawiera wyposażenie niezbędne do przeprowadzenia procedur chirurgicznych.



**Nr Kat. 4614**

**DUŻY ZESTAW STARTOWY**



Nr Kat. produktu	Opis produktu	ilość
117	Healing Abutment L4.0mm	2
4055	Hex Driver 1.25mm	1
4080	Parallel Depth Guide	3
4100	Implant Depth Probe	1
4151	Long Hex Driver 2.5mm	1
4161	Motor Mount Hex Driver 2.5/1.25mm	1
4240	Drill Extension L17.5mm	1
<b>4550</b>	<b>Coated Stright Drill BD2.0</b>	1
<b>4552</b>	<b>Coated Stright Drill BD2.8</b>	1
<b>4553</b>	<b>Coated Stright Drill BD3.0</b>	1
<b>4554</b>	<b>Coated Stright Drill BD3.2</b>	1
<b>4555</b>	<b>Coated Stright Drill BD3.65</b>	1
<b>4556</b>	<b>Coated Stright Drill BD4.1</b>	1
<b>4557</b>	<b>Coated Stright Drill BD4.5</b>	1
<b>4558</b>	<b>Coated Stright Drill BD4.8</b>	1
4572	Universal Torque Ratchet 10-45Ncm	1
46130	Advanced Surgical Organizer Box	1
4712C	Marking Drill-Sphere shape 1.5/L16	1
5061	Open Tray Transfer	2
<b>NICE 7305</b>	<b>Implant Driver Standard 2.1mm CHC</b>	1

Ważne informacje

